



List of monoclonal antibodies specific Flavivirus

No.	Isotype	Flaviviruses								Ag. Molecule
		D1	D2	D3	D4	JE	YF	WN	KUN	
YH 0001	G2a	-	-	-	-	4+	-	-	-	NS1 conformation
YH 0002	G2a	-	-	-	-	4+	-	-	-	NS1 linear
YH 0003	G1κ	4+	-	-	-	-	-	-	-	NS1 conformation
YH 0004	G1	-	3+	-	-	2+	-	-	-	NS1 monomer
YH 0005	G1	4+	2+	-	-	1+	-	-	-	NS1 dimer
YH 0006	G2 bk	-	3+	2+	-	1+	-	-	-	E conformation
YH 0007	Mκ	-	-	3+	-	-	-	-	-	NS1 conformation
YH 0008	G1κ	-	-	-	4+	-	-	-	-	NS1 conformation
YH 0009	G2aκ	-	-	-	-	3+	-	3+	3+	E conformation
YH 0010	G1	-	3+	1+	-	1+	2+	±	±	NS1 conformation
YH 0011	G2aκ	2+	3+	1+	3+	1+	-	-	-	NS1 conformation
YH 0012	G2aκ	-	1+	1+	1+	4+	-	±	-	E conformation
YH 0013	G1	-	1+	-	-	4+	-	3+	-	NS1 conformation
YH 0014	G1	-	3+	-	-	4+	-	±	-	NS1 conformation
YH 0015	G1	-	3+	-	-	2+	-	-	-	NS1 conformation
YH 0016	G2a	4+	4+	-	-	4+	-	3+	2+	NS1 conformation
YH 0017	G2a	-	-	-	-	4+	-	2+	2+	E conformation
YH 0018	G1κ	-	-	-	-	4+	2+	3+	3+	NS1 conformation
YH 0019	G1κ	-	-	-	-	4+	2+	4+	3+	NS1 conformation
YH 0020	G2aκ	-	-	-	-	-	-	4+	3+	NS1 conformation
YH 0021	G1 κ	-	4+	-	-	3+	-	-	-	E conformation
YH 0022	Mκ	1+	2+	1+	1+	2+	-	-	-	E conformation
YH 0023	G1	4+	4+	4+	4+	4+	4+	4+	4+	NS1 linear
YH 0024	G1λ	4+	3+	3+	4+	4+	4+	3+	2+	NS1 linear
YH 0025	G2aκ	3+	3+	3+	3+	3+	4+	3+	3+	E conformation
YH 0026	G1λ	4+	4+	4+	4+	4+	4+	4+	4+	E conformation
YH 0028	G2aκ	1+	-	-	-	-	-	4+	1+	E linear
YH 0029	G2aκ	-	-	-	1+	2+	1+	2+	2+	E linear
YH 0030	G2aκ	-	-	-	-	-	4+	-	-	prM linear
YH 0031	G2aκ	-	-	-	-	-	3+	-	-	prM linear
YH 0032	G3κ	-	3+	-	-	-	-	-	-	NS1 conformation
YH 0033	G1κ	3+	2+	3+	3+	3+	3+	3+	3+	NS2(A+B) linear
YH 0034	G3κ	3+	3+	3+	3+	3+	-	3+	3+	NS3 conformation
YH 0035	G2bk	2+	3+	-	±	3+	-	1+	-	E conformation
YH 0036	G1κ	2+	2+	2+	2+	2+	2+	2+	2+	E conformation





List of monoclonal antibodies specific Flavivirus

No.	Isotype	Flaviviruses								Ag. Molecule
		D1	D2	D3	D4	JE	YF	WN	KUN	
YH 0037	G1κ	2+	2+	2+	2+	2+	2+	2+	2+	E conformation
YH 0038	G1	-	-	-	-	4+	-	-	-	NS1 conformation

★ 專一辨識推薦

NO.	Isotype	Flaviviruses								Ag. Molecule
		D1	D2	D3	D4	JE	YF	WN	KUN	
YH0001	G2a	-	-	-	-	4+	-	-	-	NS1 conformation
YH0002	G2a	-	-	-	-	4+	-	-	-	NS1 linear
YH0003	G1κ	4+	-	-	-	-	-	-	-	NS1 conformation
YH0007	Mκ	-	-	3+	-	-	-	-	-	NS1 conformation
YH0008	G1κ	-	-	-	4+	-	-	-	-	NS1 conformation
YH0030	G2aκ	-	-	-	-	-	4+	-	-	prM linear
YH0032	G3κ	-	3+	-	-	-	-	-	-	NS1 conformation
YH0038	G1	-	-	-	-	4+	-	-	-	NS1 conformation

★ 通用辨識推薦

NO.	Isotype	Flaviviruses								Ag. Molecule
		D1	D2	D3	D4	JE	YF	WN	KUN	
YH0023	G1	4+	4+	4+	4+	4+	4+	4+	4+	NS1 linear
YH0026	G1λ	4+	4+	4+	4+	4+	4+	4+	4+	E conformation

★ 感謝各研究院所老客戶誠心推薦

NO.	Isotype	Flaviviruses								Ag. Molecule
		D1	D2	D3	D4	JE	YF	WN	KUN	
YH0001	G2a	-	-	-	-	4+	-	-	-	NS1 conformation
YH0003	G1κ	4+	-	-	-	-	-	-	-	NS1 conformation
YH0007	Mκ	-	-	3+	-	-	-	-	-	NS1 conformation
YH0008	G1κ	-	-	-	4+	-	-	-	-	NS1 conformation
YH0023	G1	4+	4+	4+	4+	4+	4+	4+	4+	NS1 linear
YH0025	G2aκ	3+	3+	3+	3+	3+	4+	3+	3+	E conformation
YH0026	G1λ	4+	4+	4+	4+	4+	4+	4+	4+	E conformation
YH0034	G3κ	3+	3+	3+	3+	3+	-	3+	3+	NS3 conformation

- 1、以上單株抗體皆分別經 Flaviviruses 病毒株 (D1、D2、D3、D4、JE、YF、WN、KUN) 之嚴謹力價測試，並以1+~4+ (最強) 表示強度。
- 2、曜鴻採用最高品質的 Immunodeficiency (SCID) mice 專用鼠，量產抗體。
- 3、抗體應用廣：可用於 ELISA、IHC or Western Blot。
- 4、超值大包裝：1,000 μL / per set 特價優惠中。

